

Comparative Evaluation of Filtration and Imaging Properties of Filter Membranes for Microplastic Capture and Analysis

Jared Carter¹, Erin Butler¹, Corinne Amato¹, Teagan Horan¹, Gregory Madejski² and James Roussie^{1*}


1. SiMPore Inc., West Henrietta, NY 14586
2. Parverio Inc., West Henrietta, NY 14586
* Corresponding Author: jroussie@simpore.com

Objectives

Pollution by microplastics (MPs) in drinking, fresh and ocean waters, as well as food and beverages, is a growing problem that is now well-recognized in both the popular media and scientific literature. Concerning levels of MPs have been found in food and local waters, as well as in human tissues. Several government regulatory bodies have, or are in the process of, implementing mandatory MP drinking water testing. Consequently, there are needs for greater understanding of the performance characteristics of common MP analytical methods and for standardizing methods and reporting.


Here, we report on the comparative evaluation of the filtration and imaging properties of five filter membranes capable of MP capture and analysis. This study was undertaken as part of an inter laboratory methods evaluation study coordinated by the Southern California Coastal Water Research Project. We compared track-etched polycarbonate +/- gold coating (PCTE and PCTG), polytetrafluoroethylene (PTFE) porous silicon (PS) and gold-coated microslit silicon nitride membranes (MSSN-Au). Four of the filter types were purchased with a nominal 1.0 μm cut-off, except for PCTG which was purchased with a 0.8 nominal cut-off. We further demonstrate use of MSSN filters to monitor MP entrainment along a municipal drinking water delivery route.

Filter Characteristics




Filter Type	Vendor	CAT	LOT	Composition	Pore Diameter (Cutoff, microns)	Porosity (%)	Thickness (micron)	Retail Price/Filter (USD, lowest pack size)
PTFE	Advantec	H100A025A	90227640	Polytetrafluoroethylene – hydrophilic treatment	1.0	UNK	35	\$2.67 ^a
PCTE	Sterlitech	PCT1025100	M-180227	Polycarbonate – Track Etched	1.0	16%	11	\$0.91 ^b
PCTG	Sterlitech	1270007	TPA.PC-20 01.1.8	Polycarbonate – Track Etched & 40nm Gold Coated	0.80	15%	9	\$23.84 ^c
Porous Si	Smart Membranes (Thermo)	950789-W16	UNK	Macroporous Silicon Lift-Off Membrane	1.0	~60%	220	\$30 (Thermo) ^d \$19.30 (Smart Membranes) ^e
MicroSlit SiN	SiMPore	MSSN400-1.0-3LAUJ	4882	Silicon Nitride & 120nm Gold Coated	1.0	16%	0.40	\$16 ^f

Capture and Analysis Workflow [1]




Representative Microscopy



Filter Type	Phase (Exposure [ms] /Gain)	Fluorescence (Exposure [ms] /Gain)
PTFE	350/5	350/5*
PCTE	350/5	776/9
PCTG	350/5	840/9
Porous Si	N/A	N/A
Microslit-SiN	350/5	73/22


- PTFE Background fluorescence* prevented use of higher gain to improve contrast
- PS filter data not available due to impractically slow sample processing and incompatibility with phase microscopy
- Metal coating of polymeric filters tends to improve microscopy results with considerable increase in sample processing durations and filter cost
- Solid state materials (e.g., PS, Al oxide, etc.) currently in use for IR and Raman improve fluorescent microscopy results but with substantial decrease in sample processing rate and considerable expense

Clean Water Flux





- Polymeric Membranes (PTFE, PCTE) offer similar flow rates for the same surface area
- Addition of gold to PCTE (i.e., PCTEG) greatly reduces the water flow rate (~5x)
- PS filters are impractically slow (> 2 hrs for 1 L of clean water)
- MSSN-Au filters (25 mm disc format) demonstrated the fastest filtration times (>17% faster than PTFE) using substantially less filter area (6.3 vs 70.9 mm²)

Survey of Municipal Water [1]



- We used MSSN filters to survey the MP concentration along a municipal water route (Hemlock Lake, NY).
- Based on consuming 500 mL of water a day, a consumer at a final delivery site (a drinking fountain on the University of Rochester campus) would consume 496 MP particles/day and 3–7 μg of MP upon immediate use.
- If unfiltered, the pipes directly below the drinking fountain indicate 5.0 particles/mL (2522 particles/day) and 0.5–1.2 μg of MP load.

Analytical Optionality



- Simple filtration appears to be effective in reducing the debris load, as there is a ~50% reduction in the amount of particulate at the drinking fountain immediately compared to the building's source. However, more of these particulates appear to be plastic.
- Many of the observed particles in pipes before the drinking fountain are rust and sand from the environment in which the water resides.

Conclusions

- Polymeric filters generally offered good flow rates for clean water samples at reasonable cost, but at the expense of poor overall microscopy performance, owing to their heterogeneous internal structure (visible by SEM) and their intrinsic chemical composition.
- Non-polymeric MSSN filters lacking appreciable internal structure offer higher filtration rates than polymeric filters and demonstrate microscopy performance similar to comparable solid state filters.
- Future plans include expanding the survey to other available filter types and broader samples, as well as perform on-filter sample digestion to align more closely with current and future optimized analytical workflows for microplastics analysis.

Acknowledgements

[1] - Gregory R. Madejski, S. Danial Ahmad, Jonathan Musgrave, Jonathan Flax, Joseph G. Madejski, David A. Rowley, Lisa A. DeLouise, Andrew J. Berger, Wayne H. Knox, and James L. McGrath, "Silicon Nanomembrane Filtration and Imaging for the Evaluation of Microplastic Entrainment Along a Municipal Water Delivery Route", Sustainability 2020, 12(24), 10655; <https://doi.org/10.3390/su122410655>

Conflicts of Interest: Gregory Madejski is a cofounder of Parverio Inc. James Roussie is a co-founder of SiMPore Inc.

The authors are grateful for support of NIEHS Grant No. 2-R44-ES031036-02. Material and programmatic support was partially provided by Dr. Charles Wong of the Southern California Coastal Water Research Project, Filter Augmentation Study.