

Development of a Novel 20 μm Cut-off Microporous Silicon Nitride Membrane for Separating and Analyzing Microplastic Particles in Potable Water

Teagan Horan¹, Joshua Miller¹
 Jared Carter¹, Gregory Madejski², Michael Deible³ and James Roussie^{1*}

1. SiMPore Inc., West Henrietta, NY 14586
 2. University of Rochester, NY 14620
 3. RJ Lee Group, Monroeville, PA 15146
 * Corresponding Author: jrroussie@simpore.com



Abstract

Track-etched polycarbonate (PCTE) membranes have been used as a filtration standard for microplastic (MP) capture and analysis across a variety of studies.

However, the track-etch method creates pores of varying angles that complicate on-membrane particle analysis and sometimes create merged multi-pores that allow passage of MPs > 20 μm . Additionally, the membrane's overall $\sim 3 \mu\text{m}$ thickness makes it prone to folding and wrinkling, which may complicate microscopy analysis.

To address these issues, a novel 20 μm gold-coated microporous silicon nitride (MPSN-Au) membrane was developed and compared directly against gold-coated PCTE (PCTE-Au) membranes via manual manipulation and processing time, pore characteristics, light microscopy particle analysis, and Raman/FTIR analysis.


We found that on an area-normalized basis, MPSN-Au membranes offered greater gas and water flux over PCTE-Au. The regular pore geometry of MPSN-Au membranes made particle imaging and spectral analysis more consistent and easier to discriminate between captured PS particles, when compared to PCTE-Au membranes.

Total handling and processing time for each membrane was compared, which determined that the total processing time (including filtration, automated image acquisition, and particle counting) was 161.56% faster on average for MPSN-Au than PCTE-Au membranes.


In an 8-hour workday, 85.13% more MPSN-Au membranes can be handled, processed, and imaged than PCTE-Au, which equates to 66 vs 26 total samples, respectively.

Overall, these data demonstrate the utility of MPSN-Au membranes and suggest they can significantly improve testing time-related efficiency in all aspects of normal use-case situations as compared to PCTE-Au membranes.


Water and Gas Flux




Manipulation Processes



Membrane Characteristics




PORE SIZE COMPARISON - MPSN-Au vs PCTE-Au



Filter type	Composition	Pore Diameter	Porosity	Thickness (μm)	Gold coating thickness (nm)	Vendor	CAT
MPSN-Au	Gold-coated silicon nitride	18.71 μm	27.7%	0.40	120	SiMPore	5065
PCTE-Au	Gold-coated polycarbonate track-etched	23.60 μm	12%	3	20	SPI Supplies	E20025-MB

Microscopy Particle Counting




Algorithmic Particle Analysis of Light Microscopy Images

Imaging and particle analysis system - Keyence vhx-7000

- MPSN-Au - appeared to have larger overall particle sizes captured on the membrane (18.7 - 426.5 μm range).
- Resultant particle counts - 1810
- PCTE-Au - Smaller particle sizes counted, many of which appear to pores mistaken for particles based on diameter sizes returned (2.8 - 166 μm range)
- Resultant particle counts - 9382

Raman and FTIR Analysis

Raman Analysis of Blanks




Polystyrene sourced from from HPU's Polymer Kit 1.0, (J. Lynch, NIST). Data uploaded to OpenSpecy and coefficients of correlation analyzed:
 MPSN-Au - 0.77, 0.73, 0.73
 PCTE-Au - 0.75, 0.81, 0.79


Raman
 Horbia LabRAM HR with LabSpec v6 was used for all analyses.
 • 632 nm Laser, $\sim 10 \text{ mW}$ power. Spot size 1028 nm with 50x long working distance (LWD) objective.
 • 600 gr/mm grating gives $\sim 2.7 \text{ cm}^{-1}$ resolution.

FTIR
 Bruker Tensor 27 with Hyperion microscope attachment used for all analyses.
 • Samples analyzed in reflection mode.
 • 32 accumulations used for samples.

Raman Analysis of PS Particles



FTIR Analysis of Substrate



Conclusions

- MPSN-Au membranes offer a faster and more consistent method of manual manipulation than PCTE-Au.
 - When accounting for experimental filtration times, automated acquisition, and manual manipulation time, 85.13% more MPSN-Au membranes can be processed than PCTE-Au (66 vs 26, respectively, in an 8-hour workday)
- MPSN-Au had 137.34% and 151.54% greater surface area-normalized water and gas flux rates, respectively, versus PCTE-Au, but an experimentally 120.25% slower filtration rate due to the smaller surface area.
- MPSN-Au offers a highly consistent pore geometry, which improves particle counting methods by eye and by algorithm. PCTE-Au's inconsistent pore geometry and angle may lead automated counting systems to count pore edges as particles and misrepresent particle counts by a significant margin. However, the concentration factor of the MPSN-Au membranes may result in multilayer formation with the same filtered volume due to the smaller surface area.
- PCTE-Au has higher coefficient of correlation with returned Raman spectra by 5.24% on average.
- PCTE-Au is not suitable for FTIR analysis at this gold coating thickness - the PCTE substrate can be seen in the resultant spectra, while MPSN-Au returns a consistent and predictable background spectra.

Acknowledgements

- Funding Support from NIH Grant No. NIEHS 2-R44-ES031036-02, awarded to SiMPore Inc.
- SEM collection conducted at Integrated Nanosystems Center (URnano) at the University of Rochester.
- Guidance regarding Raman spectroscopy from Dr. Samantha Romanick and Dr. Andrew Berger, University of Rochester.
- Microfabrication was carried out at the Semiconductor and Microsystems Fabrication Laboratory (SMFL) at the Rochester Institute of Technology.

