


**SiMPore**  
Precision Membrane Technologies™

# μSiM-EV: Silicon Membrane-Enabled Microfluidics for Simplified Single Extracellular Vesicle Visualization

SiMPore Inc., West Henrietta, NY 14586 • www.simpore.com • info@simpore.com




© SiMPore Inc. ALL RIGHTS RESERVED.

## Overview


- ◆ The μSiM-EV enables imaged-based analysis of individual EVs.
- ◆ Works with EVs prepared by ultracentrifugation, size-exclusion chromatography, or filtration from a variety of biofluids.
- ◆ Requires simple pipette-driven loading and conventional epifluorescent microscopy, with no specialized equipment needed.

## μSiM-EV Concept




**Figure 1: The EV Visualization Concept.** A) Representative nanomembrane with EV-sized pores that capture individual fluorescent affinity-tagged EVs; fluorescent color indicates EVs positively labeled with affinity tags. B) Microfluidic device enabled by a silicon membrane for EV visualization ("The μSiM-EV"), with two injection ports, center well, and bottom channel with imaging-compatible layer.

## Fluorescent Bead Visualization




**Figure 4: Imaging via Fluorescence and Electron Microscopy.** A) Fluorescent polystyrene beads (45 nm Nile Red on the nanoporous membrane, 830 nm Jade Green on the microporous membrane) were injected into the μSiM-EV and imaged via fluorescence microscopy. B) The same μSiM-EV units as in A) were prepared for scanning electron microscopy (SEM) by removal of the bottom layer and depositing ~7 nm Au, then imaged via SEM.

## Visualize via Antibody Staining



**Figure 2: Immunolabeling Analysis by μSiM-EV of EVs Prepared by Dual Mode Chromatography.** Human plasma (0.5 mL) was applied to a column comprising 10 mL of Sepharose CL-4B and 2 mL FractoGel-sulfate A) CFSE (cyan) and anti-CD63 (magenta) fluorescent micrographs are shown via 1% representative fields-of-view from the μSiM-EV's membrane. B) Particle counts are rendered from ImageJ ComDet using 5 px max separation, 3 px/500 nm max particle size, and 3 intensity threshold. Images in A obtained from a Zeiss Z1 epifluorescent microscope, 40X objective.


## Count EVs




**Figure 3: Count EVs prepared by Dual Mode chromatography on the μSiM-EV.** Prepared EVs, as in Figure 2, were labeled with the amine-reactive pan-label CFSE. A) CFSE fluorescent micrograph B) Particle counts of CFSE-positive particles counted by ImageJ ComDet using 5 px max separation, 3 px/500 nm max particle size, and 3 intensity threshold.

## EV Visualization Workflow

### Injecting the Device



### Imaging the Device



Visit us at Booth #113!



SCAN ME!  
To save 10% on your first order