

# Flow Cell Development Kit

Tangential and transmembrane flow cell assembly

**Product Nos. GASKET-CS**

**Customer Service**

[info@simpore.com](mailto:info@simpore.com)

phone: 585-214-0585

**Order Information**

[www.SiMPore.com](http://www.SiMPore.com)

[sales@SiMPore.com](mailto:sales@SiMPore.com)

phone: 888-323-NANO

fax: 888-249-2935



[www.SiMPore.com](http://www.SiMPore.com)

150 Lucius Gordon Drive | Suite 110 | West Henrietta NY 14586

## INSTRUCTIONS FOR USE

### Provided Kit Materials

- Pre-Cut PDMS Gasket Set (suitable for 10 devices)
- 18Ga 1/2" Stainless Steel Luer Lock straight Needles (2)
- 18Ga 1/2" Stainless Steel Luer Lock ninety degree Needles (2)
- #1 22 x 22 mm cover slips (10)



Figure 1 (A) Gasket bag components and (B) internal component bags

### Required Equipment & Supplies

#### *Equipment*

- UV Ozone generator or Oxygen Plasma Oven
- Hot Plate or gravity oven

#### *Supplies*

- Tweezers suitable for membrane manipulation (SiMPore K6TWZR, or equivalent)
- User-supplied tubing or fluidic interface(s)
- Petri dish, glass slide, or other clean **substrate** on which to assemble devices
- Alcohol wipes or KimWipe saturated with an alcohol solution
- Absolute Alcohol (such as IPA or Ethanol)
- Wetting Media (lab prepared)
- Gloves and other suitable PPE

## Gasket Stack Layout



**Figure 2.** Gasket stack order of assembly

## INSTRUCTIONS

Please read all **instructions** before proceeding. Refer to the gasket stack diagram above which details the order of assembly. If unfamiliar with assembling microfluidic devices from PDMS, refer to the Tips & Troubleshooting section for more guidance.


### **Device Assembly**

1. Working on a clean, particle-free work surface, remove the PDMS gaskets and cover slips from the supplied bags.
2. Clean the glass cover slip with an alcohol wipe, allow to air dry, and place on the **substrate**.
3. Using tweezers transfer the lower channel layer (1) to the substrate and transfer to a UV ozone or O<sub>2</sub> Plasma source, then surface activate as per tool specifications.

**NOTE:** Reference Figure 1 (above) for gasket layer number and order of assembly.

4. Carefully transfer the now surface activated Lower Channel Layer (1) to the cover slip, aligning the two to ensure an even fit.
5. Repeat the above process (steps 3 - 4) for adhering the Membrane Interface Layer (2) onto the Lower Channel Layer (1).
6. Continue the layer addition process as above to secure the Membrane Housing Layer (3) to the Interface Layer (2)
7. Retrieve a 5.4 mm x 5.4 mm SiMPore membrane chip and surface activate the lower assembled layers and chip as before

**NOTE:** When manipulating the membrane chips, ensure the membrane-side of the chip remains **up** at all times to prevent membrane damage. Reference Figure 3 below.



**Figure 3.** Membrane chip organization. SiMPore membranes are deposited across the surface of a silicon wafer support material. Access to the membrane is achieved by etching trenches through the silicon wafer support, yielding the cross sectional geometry shown above.

8. Install the now surface activated membrane chip into the lower channel assembly such that the **membrane** side is **down** and the long axis of the **trench** are parallel to the long axis of the gasket layers.
 

**NOTE:** Incorrect placement or orientation of the membrane chip will prevent proper device function, and likely will result in membrane failure
9. Retrieve and surface activate the upper Membrane Interface Layer (4) and surface activate with the membrane assembly, then assemble as shown in Figure 2.
10. Surface activate the Upper Channel Layer (5) and surface activate with the membrane assembly, then assemble as shown in Figure 2.
11. Finally, retrieve an Interface Layer (6) and surface activate with the membrane assembly, then assemble as shown in Figure 2.
 

**NOTE:** The interface layer is 800 microns in total thickness and must be well-aligned to ensure an open fluid path through the device.
12. Once fully assembled, cure the device for at least 60 minutes at 80°C

### **Using the Device**

1. After curing, use the provided 18Ga needles to develop fluid interfaces to the device.
  - Slowly with even pressure insert the needle through the port within the interface layer. Some resistance is expected.
  - Ensure proper placement by visual examination of the device from the side and top such that the steel needle does not embed into the bottom of the channel layer.
  - Repeat the above process for the remaining three needles.
2. Once ported, connect tubing assemblies (user supplied) to the Luer lock ports and prime using a wetting solution (particle free (0.2  $\mu\text{m}$  filtered) dilute IPA or other hygroscopic solution) by slow injection from a syringe. If resistance is

www.SiMPore.com

150 Lucius Gordon Drive | Suite 110 | West Henrietta NY 14586

experienced during priming, slowly withdraw the needle(s) to enable fluid flow through each flow channel.

**NOTE:** If significant transmembrane flow is observed, discontinue flow and examine the outlet channel and port for closures.

3. Continue priming with wetting solution until the entire flow path is purged of air.

**NOTE:** Small bubbles encapsulated within the device may be dislodged by gentle tapping. Alternatively, slowly withdrawing the solution front through the device will cause small bubbles to pop against the retreating solution-air interface.


4. Once fully wetted, exchange the wetting solution with the desired media (experiment specific) with sufficient volume to ensure displacement of all residual wetting solution
5. As a general rule, avoid experimental parameters that exceed the following to ensure device and membrane integrity is not compromised

Parameter	Value
Flow Rate	<5 mL/min
Membrane $\Delta P$	< 68 kPa
Device $\Delta P$	< 55 kPa
Media Composition	Aqueous solutions, solvents compatible with PDMS. De-gassed and particle-free

## TIPS AND TROUBLESHOOTING

- **Inspection:** Prior to assembling device layers, carefully inspect the gaskets for debris, which may include PDMS cuttings or plugs (particularly in fluid ports or other small feature).
- **Cleanliness:** It is of critical importance to maintain the cleanliness of the PDMS surface to ensure proper bonding between gasket layers. If dust, debris, or other visible debris is present, attempt to remove via tweezers or compressed nitrogen gas. If accidental contact is made with an ungloved hand the PDMS gasket must be thoroughly cleaned with IPA or other solvent to remove any surface oil.
- **Surface Activation:** Surface activation of PDMS should be performed via UV Ozone or PDMS. As tool settings and operating parameters vary, development of an optimized recipe for complete surface activation is left to the user. Once surface activated the PDMS should be bonded as quickly as possible, within five minutes for optimal bonding results.
- **Proper Sealing:** Good bonding between the two PDMS layers is critical to ensure device performance. Debris or bubbles between the layers may cause device failure. Once bonded, no movement between the layers is advised (e.g. replacement during alignment), as the bond strength is greatly diminished with repeated contact. If misplacement occurs, repeat the surface activation before bonding again.

- **Curing:** Once fully assemble, the device must be heat-cured before use to improve the integrity of the bonded layers. Bake the assembled device on a hot plate or more ideally a gravity oven, for at least 60 minutes at 80°C. Bake times will vary depending on the tool used.
- **Device Priming:** Once cured and the needles inserted, the device layer and membrane may be initially hydrophobic due to off-gassing during the cure. A hygroscopic solution, ideally dilute alcohol or glycerol in MilliQ water may be used to fully wet the channels and membrane before use. Rinse the device fully with the desired contact media before experimentation.
- **Using the Device:** Adhere to the recommended values provided above in *Using the Device*, Step 5. Caking or fouling may occur if particle or biofouling agent density is too high leading to device pressure increase and membrane or device failure. Device service life has not been fully tested and varies depending on assembly quality and use case(s), but generally is suitable for at least 72 hours of continuous use with diligent care. Refer to Figure 4 below for a flow path diagram and picture of a represented device after assembly.



**Figure 4 (A)** Assembled gasket geometry stack and direction of flow for the upper (red) and lower (blue) channels. **(B)** A representative image of a fully assembled flow cell containing a membrane chip

## PHYSICAL PROPERTIES

Property	Value
Membrane Composition:	Silicon Nitride (SiN)
Membrane Surface Area	1.4 mm <sup>2</sup> (single slot), 6.3 mm <sup>2</sup> (three slot)
Minimum Working Length	750 μm
Membrane Thickness:	100 nm (NPN) / 400 nm (SiN)
Pore Size Cut-Off:	varies
Porosity:	varies
Surface Charge:	Neutral to slightly negative
Wetted Components Material(s):	PDMS, Si, SiN, SS

## SPECIFICATIONS

### Sterilization

Flow Cell devices may be sterilized via Ethylene Oxide, Hydrogen Peroxide Gas, 70% IPA immersion, or gamma, UV, and E-Beam irradiation after removing from the plastic bag. **Do Not** Sterilize via steam autoclave as device damage may occur.

### Chemical Stability

Devices are incompatible with strong bases, as membrane degradation may occur. Avoid long exposure durations to solvents, as this may compromise device integrity. Do not apply solvents, acids, bases, organics, or other compounds that are incompatible with the stated device components.

### Device Storage

Store in a clean and dry environment. Prolonged exposure to UV can degrade the PDMS bonding, leading to poor layer adhesion and shortened service life.

---

## STANDARD WARRANTY

**SiMPore Inc.** ("SiMPore") warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. **SIMPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications, and descriptions of SiMPore products appearing in SiMPore's published catalogues and product literature may not be altered except by express written agreement signed by an officer of SiMPore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, SiMPore's sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies SiMPore promptly of any such breach. If after exercising reasonable efforts, SiMPore is unable to repair or replace the product or part, then SiMPore shall refund to the customer all monies paid for such applicable product or part. **SIMPORE SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL, OR ANY OTHER INDIRECT DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY CUSTOMER FROM THE USE OF ITS PRODUCTS.**

## TRADEMARKS AND PATENTS

This product is covered by US Patents # 9,789,239 & #10,391,219

Version 1226