

How to Choose a TEM Window Grid

Application Note: Guidelines for selecting a TEM Window Grid

Slots or squares?

Our TEM windows are available in an assortment of geometries including standard two slot or nine square formations. Due to the geometry, the squares will be slightly more robust and more likely to survive rigorous sample preparation. The slots allow larger uninterrupted viewing windows

100 or 200 micron thick frames?

All of our grids are designed to fit into standard 3 mm TEM sample holders. Our 200 micron thick frames fit into most standard holders. Additionally, we produce 100 micron thick frames that fit holders designed for high tilt or thin carbon grid applications.

Membrane thickness and type?

The choice of membrane thickness and type is dependent on many variables, including: Thinner membranes (5–10 nm) are more suitable for higher resolution imaging, while thicker membranes (15–50 nm) are better for demanding sample preparation procedures.

Consider using a thicker membrane for method development and moving to thinner membranes for final imaging. The nature of the specimen preparation may determine the choice of membrane. The thermal and chemical stability, surface hydrophobicity, and plasma cleanability differ among the various membranes and should be taken into consideration for each application.

	Amorphous Silicon	Porous Nanocrystalline Silicon	Silicon Dioxide	Silicon Nitride	Standard Carbon	Ultrathin Carbon
Actual Thickness (nm)	5, 9 & 15	15	20 & 40	10, 20 & 50	20–50	~10
Image Quality	Excellent	Good	OK	Good	OK	Good
Plasma Cleanable	Yes	Yes	Yes	Yes	No	No
Elemental Analysis Background	Si Only	Si Only	Si, O	Si, N	C, H	C, H
Thermal Stability	~600 C	>1000C	>1000C	>1000C	~400 C	~400 C
Chemical Stability	Avoid Strong Bases	Avoid Strong Bases	Good	Excellent	Good	Good
Tolerates High Beam Currents	Excellent	Excellent	OK	OK	Excellent	Excellent
Potential Contamination Source	None	None	None	None	Carbon	Carbon
Open Nanoscale Pores	No	Yes	No	No	No	No
Background	Featureless	Nanocrystalline	Featureless	Featureless	Featureless	Featureless